

ABLATION THERAPY

Ablation Catheter Curve Reference Charts.	A.1
--	-----

Standard Ablation Catheters

Livewire™ TC Ablation Catheters.	A.13
Therapy™ Ablation Catheters.	A.16
Therapy ComfortGrip™ Bi-directional Ablation Catheter	A.22

Irrigation Ablation Catheters

FlexAbility™ Ablation Catheter, Sensor Enabled™	A.23
FlexAbility™ Ablation Catheters	A.24
TactiCath™ Contact Force Ablation Catheter, Sensor Enabled™	A.25
TactiCath™ Quartz Contact Force Ablation Catheter	A.26
Therapy™ Cool Path™ Ablation Catheters	A.27
Therapy™ Cool Path™ Duo Ablation Catheters	A.29
Therapy™ Cool Flex™ Ablation Catheter	A.30
Contact™ Therapy™ Cool Path™ Duo Ablation Catheter	A.31

Hardware




Cool Point™ Irrigation Pump	A.32
Ampere™ RF Ablation Generator	A.33
TactiSys™ Quartz Equipment	A.34

Connecting Cables

Connecting Cables for Livewire™ TC Ablation Catheters	A.35
Connecting Cables for Therapy™ and Contact Enabled™ Ablation Catheters	A.36

Livewire™ TC Ablation Catheter


(4 mm bi-directional)
Curve Reference



Livewire™ TC Ablation Catheter

(4 & 8 mm uni-directional)



Curve Reference



Livewire™ TC Ablation Catheter

(4 & 8 mm uni-directional)

Curve Reference




Therapy™ Ablation Catheter

(4 & 8 mm uni-directional) and


Therapy™ Dual-8™ Ablation Catheter

(8 mm uni-directional) Curve Reference


SMALL CURVE



MEDIUM CURVE




LARGE CURVE




TactiCath™ Quartz Ablation Catheter

Curve Reference

UNI-DIRECTIONAL 65



UNI-DIRECTIONAL 75



FlexAbility™ Ablation Catheter

Bi-Directional Curve Reference


FlexAbility™ Irrigated Ablation Catheter,

Sensor Enabled™


TactiCath™ Contact Force Ablation Catheter,

Sensor Enabled™ Curve Reference


BI-DIRECTIONAL CURVE D-D




BI-DIRECTIONAL CURVE D-F




BI-DIRECTIONAL CURVE F-F



BI-DIRECTIONAL CURVE J-J



BI-DIRECTIONAL CURVE F-J



FlexAbility™ Ablation Catheter

Uni-Directional Curve Reference


FlexAbility™ Irrigated Ablation Catheter,

Sensor Enabled™


TactiCath™ Contact Force Ablation Catheter,

Sensor Enabled™ Curve Reference


UNI-DIRECTIONAL CURVE D



UNI-DIRECTIONAL CURVE F



UNI-DIRECTIONAL CURVE J






Therapy™ Ablation Catheter

(4 & 8 mm) and

Therapy™ Dual-8™ Ablation Catheter

(8 mm uni-directional)

Curve Reference




Therapy™ Ablation Catheter


(4 & 8 mm bi-directional)

Curve Reference


SMALL CURVE




MEDIUM CURVE



LARGE CURVE








X-LARGE CURVE




Therapy™ Cool Path™ Irrigated Ablation Catheter

(2 & 4 mm uni-directional)
Curve Reference







Therapy™ Cool Path™ Duo, Therapy™ Cool Flex™ and Contact™ Therapy™ Cool Path™ Duo Ablation Catheters

4 mm Tip
Curve Reference



Therapy™ Cool Path™ Duo Ablation Catheters (Bi-directional)

4 mm Tip
Curve Reference



Livewire™ TC Ablation Catheters

4 mm Tip Universal
Temperature Monitoring
Quadripolar
7 F

Product Highlights

- Bi-Directional tip deflection
- Thermocouple and Thermistor temperature monitoring capability
- Rotary dial designed for fine tip movements
- Ergonomic handling

Ordering Information


REORDER NUMBER	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	DISTAL REACH	USABLE LENGTH (CM)
402275	2-5-2	4	1	Med Sweep	115
402277	2-5-2	4	1	Med Curl	115
402278	2-5-2	4	1	Lrg Curl	115

Required Catheter Extension Cables – Page A.35

SJM: 402543, 402544

Stockert: 402505

EPT: 402542



Livewire™ TC Ablation Catheters

4 & 8 mm Tip Universal
Temperature Monitoring
Quadripolar
7 F

Product Highlights

- Thermocouple and Thermistor temperature monitoring capability
- Rotary dial designed for fine tip movements
- Extra lateral stability (XLS) incorporated at catheter tip
- Ergonomic handling
- Uni-directional tip deflection

Ordering Information

7 F Quadripolar (XLS) catheter (1 unit per box)


REORDER NUMBER	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	DISTAL REACH	USABLE LENGTH (CM)
402112	2-5-2	4	1	Med Sweep	115
402113	2-5-2	4	1	Lrg Sweep	115
402114	2-5-2	4	1	Med Curl	115
402115	2-5-2	4	1	Lrg Curl	115
402193	2-5-2	8	1	Med Sweep	115
402194	2-5-2	8	1	Lrg Sweep	115
402195	2-5-2	8	1	Med Curl	115
402196	2-5-2	8	1	Lrg Curl	115

Required Catheter Extension Cables – Page A.35

SJM: 402543, 402544

Stockert: 402505

EPT: 402542



Livewire™ TC Uni-directional Ablation Catheters

4 & 8 mm Tip Thermocouple
Quadripolar
7 F

Product Highlights

- Thermocouple temperature monitoring capability
- Rotary dial designed for fine tip movements
- Ergonomic handling
- Uni-directional tip deflection


Ordering Information

REORDER NUMBER	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	DISTAL REACH	USABLE LENGTH (CM)
402120	2-5-2	4	Sml Sweep	115
402121	2-5-2	4	Med Sweep	115
402122	2-5-2	4	Lrg Sweep	115
402124	2-5-2	4	Med Curl	115
402125	2-5-2	4	Large Curl	115

Required Catheter Extension Cables – Page A.35

SJM: 402543, 402544

Stockert: 402505



Therapy™ Ablation Catheters

4 mm Tip Thermocouple
5 F

Product Highlights

- Uni-directional steering
- Thermocouple temperature measurement
- Push/pull handle

Ordering Information

5 F Quadripolar catheter (1 unit per box)


REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE
IBI-83351	1304-5-25-M-TE4BE1 (SOFT)	2-5-2	4	1	Medium
IBI-83403	1304-5-2-S-TE4BE1	2	4	1	Small

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641

Stockert: IBI-85713, IBI-85709

Medtronic: IBI-85711, IBI-85708



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermocouple
Quadripolar
7 F

Product Highlights

- Uni-directional steering
- Thermocouple temperature measurement
- Push/pull handle

Ordering Information

7 F ablation catheter (1 unit per box)


REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83302	1304-7-25-M-TE8	2-5-2	8	2	Medium	110
IBI-83306	1304-7-25-L-TE8	2-5-2	8	2	Large	110
IBI-83309	1304-7-25-L-TE8A	2-5-2	8	2	Large	110
IBI-83311	1304-7-25-XL-TE8A	2-5-2	8	2	X-Large	110
IBI-83312	1304-7-25-XL-TE8	2-5-2	8	2	X-Large	110
IBI-83404	1304-7-2-M	2	4	2	Medium	110
IBI-83405	1304-7-25-M	2-5-2	4	2	Medium	110
IBI-83408	1304-7-25-L	2-5-2	4	2	Large	110
IBI-83411	1304-7-25-E	2-5-2	4	2	Extended Reach	110
IBI-83417	1304-7-25-XL	2-5-2	4	2	X-Large	110
IBI-83432	1304-7-25-S	2-5-2	4	2	Small	110

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641

Stockert: IBI-85713, IBI-85709

Medtronic: IBI-85711, IBI-85708



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermistor

Quadripolar

7 F

Product Highlights

- Uni-directional steering
- Thermistor temperature measurement
- Compatible with EPT-1000 XP generator
- Push/pull handle


Ordering Information

7 F ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83453	1304-7-25-S-TH	2-5-2	4	2	Small	110
IBI-83456	1304-7-25-M-TH	2-5-2	4	2	Medium	110
IBI-83459	1304-7-25-L-TH	2-5-2	4	2	Large	110
IBI-83477	1304-7-25-L-TE8-TH	2-5-2	8	2	Large	110
IBI-83481	1304-7-25-XL-TE8-TH	2-5-2	8	2	X-Large	110

Required Catheter Connecting Cable – Page A.36

EPT: IBI-85763



Therapy™ Dual-8™ Ablation Catheters

8 mm Tip Dual Thermocouple
Quadripolar
7 F

Product Highlights

- Uni-directional steering
- 8 mm tip electrode with two thermocouples
- High power option up to 100 watts
- Push/pull handle

Ordering Information


7 F ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83421	1304-7-2-M-TE8TC2	2	8	2	Medium	110
IBI-83422	1304-7-25-M-TE8TC2	2-5-2	8	2	Medium	110
IBI-83424	1304-7-2-L-TE8TC2	2	8	2	Large	110
IBI-83425	1304-7-25-L-TE8TC2	2-5-2	8	2	Large	110

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641

Stockert: IBI-85739



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermocouple Triflex
 Quadripolar
 7 F

Product Highlights

- Uni-directional steering
- 4 mm tip with one thermocouple; 8 mm tip with two thermocouples
- Push/pull handle
- Three shaft sections with gradually increasing softness (braided proximal shaft, semi-soft middle shaft, soft deflectable tip section)

Ordering Information

7 F ablation catheter (1 unit per box)


REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83510	1304-7-25-M-TF	2-5-2	4	2	Medium	110
IBI-83513	1304-7-25-S-TF	2-5-2	4	2	Small	110
IBI-83516	1304-7-25-L-TF	2-5-2	4	2	Large	110
IBI-83542	1304-7-25-L-TE8TC2-TF	2-5-2	8	2	Large	110
IBI-83543	1304-7-25-XL-TE8TC2-TF	2-5-2	8	2	X-Large	110
IBI-83545	1304-7-25-M-TE8TC2-TF	2-5-2	8	2	Medium	110
IBI-83552	1304-7-25-M-TH-TF(90)	2-5-2	4	2	Medium	90

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641

Stockert: IBI-85713, IBI-85709

Medtronic: IBI-85711, IBI-85708



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermocouple Triflex
Quadripolar
7 F

Product Highlights

- Bi-directional steering
- Thermocouple temperature measurement
- Push/pull handle
- Three shaft sections with gradually increasing softness
(braided proximal shaft, semi-soft middle shaft, soft
deflectable tip section)

Ordering Information

7 F ablation catheter (1 unit per box)


REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE
IBI-83701	1304-7-25-M-BD	2-5-2	4	2	Medium
IBI-83702	1304-7-25-L-BD	2-5-2	4	2	Large
IBI-83703	1304-7-25-S-BD	2-5-2	4	2	Small
IBI-83704	1304-7-25-L-TE8-BD	2-5-2	8	2	Large
IBI-83707	1304-7-25-M-TE8-BD	2-5-2	8	2	Medium
IBI-83708	1304-7-25-XL-TE8-BD	2-5-2	8	2	X-Large

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641

Stockert: IBI-85713, IBI-85709

Medtronic: IBI-85711, IBI-85708



Therapy ComfortGrip™ Bi-directional Ablation Catheter

Product Highlights

- Ergonomic design for superb useability and efficiency
- Bi-directional deflection for improved reach in challenging anatomies
- Balanced ergonomic handle
- Intuitive steering mechanism
- Nonslip grip
- Patented steering mechanism offers smooth, precise deflection
- Braided catheter shaft provides superb torque response and stable catheter position
- Universal Temperature Monitoring, with thermocouple and thermistor sensors in every catheter
- Wide range of sweep and curl tip deflection styles
- Broad range of extension cables connect to most RF generators — simply plug in and go

Ordering Information

REORDER NUMBER	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	DISTAL REACH
A-TCG-402818	2-5-2	4	Small Sweep
A-TCG-402819	2-5-2	4	Medium Sweep
A-TCG-402820	2-5-2	4	Large Sweep
A-TCG-402821	2-5-2	4	Small Curl
A-TCG-402822	2-5-2	4	Medium Curl
A-TCG-402823	2-5-2	4	Large Curl
A-TCG-402825	2-5-2	8	Medium Sweep
A-TCG-402826	2-5-2	8	Large Sweep
A-TCG-402828	2-5-2	8	Medium Curl
A-TCG-402829	2-5-2	8	Large Curl

FlexAbility™ Ablation Catheter,

Sensor Enabled™

Product Highlights

- Combines the unique flexible tip with a handle and shaft designed to accurately navigate within the heart^{1,2}
- D, F and J curve options
- Symmetric or asymmetric options
- Ergonomic uni- and bi-directional handles are designed for consistent performance and fatigue reducing enhancements.*³
- Integrates with EnSite Precision™ cardiac mapping system and MediGuide™ Technology

Ordering Information

8 F irrigated ablation catheter

REORDER NUMBER	DESCRIPTION	CURVE	FRENCH SIZE	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	LENGTH (CM)
A-FASE-DD	Sensor Enabled™ Bi-D Irrigated Ablation Catheter	D-D	8	1-4-1	4	115
A-FASE-DF	Sensor Enabled™ Bi-D Irrigated Ablation Catheter	D-F	8	1-4-1	4	115
A-FASE-FF	Sensor Enabled™ Bi-D Irrigated Ablation Catheter	F-F	8	1-4-1	4	115
A-FASE-FJ	Sensor Enabled™ Bi-D Irrigated Ablation Catheter	F-J	8	1-4-1	4	115
A-FASE-JJ	Sensor Enabled™ Bi-D Irrigated Ablation Catheter	J-J	8	1-4-1	4	115
A-FASE-D	Sensor Enabled™ Uni-D Irrigated Ablation Catheter	D	8	1-4-1	4	115
A-FASE-F	Sensor Enabled™ Uni-D Irrigated Ablation Catheter	F	8	1-4-1	4	115
A-FASE-J	Sensor Enabled™ Uni-D Irrigated Ablation Catheter	J	8	1-4-1	4	115

Required Catheter Connecting Cable: Model A-FASE-CBL4



*Compared to Therapy™ Cool Flex™ ablation catheter

1. Abbott. Data on File. Report 90210281.

2. Abbott. Data on File. Report 90115246.

3. Abbott. Data on File. Report 90213554

The FlexAbility™ Ablation Catheter, Sensor Enabled™ is intended for use with the compatible irrigation pump and a compatible RF cardiac ablation generator. The catheter is intended for creating focal lesions during cardiac ablation procedures (mapping, stimulation, and ablation) for the treatment of arrhythmias. Sensor Enabled™ Ablation Catheters are used with compatible magnetic tracking systems to enable catheter positioning and navigation. Epicardial ablation should be limited to appropriately selected patients with ventricular tachycardia.

FlexAbility™

Ablation Catheter

Product Highlights


- Combines the unique flexible tip on the market with a handle and shaft designed to accurately navigate within the heart^{1,2}
- D, F and J curve options
- Symmetric or asymmetric options

Ordering Information

8 F irrigated ablation catheter

REORDER NUMBER	DESCRIPTION	CURVE	FRENCH SIZE	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	LENGTH (CM)
A701124	Bi-D Irrigated Ablation Catheter (4-hole)	D-D	8	1-4-1	4	115
A701125	Bi-D Irrigated Ablation Catheter (4-hole)	D-F	8	1-4-1	4	115
A701127	Bi-D Irrigated Ablation Catheter (4-hole)	F-F	8	1-4-1	4	115
A701128	Bi-D Irrigated Ablation Catheter (4-hole)	F-J	8	1-4-1	4	115
A701129	Bi-D Irrigated Ablation Catheter (4-hole)	J-J	8	1-4-1	4	115
A701157	Uni-D Irrigated Ablation Catheter (4-hole)	D	8	1-4-1	4	115
A701158	Uni-D Irrigated Ablation Catheter (4-hole)	F	8	1-4-1	4	115
A701159	Uni-D Irrigated Ablation Catheter (4-hole)	J	8	1-4-1	4	115

Required Catheter Connecting Cable: Model 1641



1. Abbott. Data on File. Report 90210281.

2. Abbott. Data on File. Report 90115246.

The FlexAbility™ Ablation Catheter is intended for use with the compatible irrigation pump and a compatible RF cardiac ablation generator. The catheter is intended for creating focal lesions during cardiac ablation procedures (mapping, stimulation, and ablation) for the treatment of arrhythmias. Epicardial ablation should be limited to appropriately selected patients with ventricular tachycardia.

TactiCath™ Contact Force

Ablation Catheter,
Sensor Enabled™

Product Highlights

- Fully integrated with the EnSite Precision™ Cardiac Mapping System and MediGuide™ Technology
- Gain additional reach by using a bi-directional handle with both symmetric and asymmetric curves, or select a uni-directional option
- Fast, easy setup and enhanced workflow is obtained through an intuitive display of CF data
- CF parameters:
 - FTI™ – Force Time Integral parameter measured in real-time during catheter ablation procedures
 - LSI™ – Lesion Index combining contact force, radiofrequency (RF) application duration, and RF current
- Fiber optic contact force technology that provides superior accuracy and reliable readings
- Requires TactiSys™ Quartz equipment and EnSite™ Contact Force Module software

Ordering Information

Open-irrigated 8 F ablation catheter

REORDER NUMBER	DESCRIPTION	CURVE	FRENCH SIZE	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	LENGTH (CM)
A-TCSE-DF	Bi-D Irrigated Ablation Catheter	D-F	8	2-2-2	3.5	115
A-TCSE-FJ	Bi-D Irrigated Ablation Catheter	F-J	8	2-2-2	3.5	115
A-TCSE-D	Uni-D Irrigated Ablation Catheter	D	8	2-2-2	3.5	115
A-TCSE-F	Uni-D Irrigated Ablation Catheter	F	8	2-2-2	3.5	115
A-TCSE-J	Uni-D Irrigated Ablation Catheter	J	8	2-2-2	3.5	115



TactiCath™ Quartz

Contact Force Ablation Catheter

Product Highlights

- Provides fast, real-time contact force (CF) measurement based on light interferometry technology
- Real-time feedback: Fast 50 Hz (20ms) sample rate streamlines workflow
- Fast, easy setup and enhanced workflow is obtained through an intuitive display of CF data
- CF parameters:
 - FTI™ – Force Time Integral parameter measured in real-time during catheter ablation procedures
 - LSI™ – Lesion Index combining contact force, radiofrequency (RF) application duration, and RF current
- Force sensor at the distal catheter tip to ensure accuracy and stability of measurement
- Fiber optic technology for optimal sensitivity and stability during CF measurement
- 75 or 65 mm steerable curve
- Requires TactiSys™ Quartz equipment and EnSite™ Contact Force Module software

Ordering Information

Open-irrigated 7 F ablation catheter with 75 or 65 mm steerable curve

REORDER NUMBER	DESCRIPTION	CURVE	FRENCH SIZE	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	LENGTH (CM)
PN-004 075	Irrigated Ablation Catheter (6 hole)	75 mm	7	2-5-2	3.5	115
PN-004 065	Irrigated Ablation Catheter (6 hole)	65 mm	7	2-5-2	3.5	115



Therapy™ Cool Path™ Ablation Catheter

Irrigated 2 & 4 mm Tip

Quadripolar

7 F

Product Highlights

- Six equidistant ports on the distal tip
- Uni-directional steering
- Push/pull handle


Ordering Information

7 F irrigated ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83352	1304-CP-7-25-M	2-5-2	4	2	Medium	110
IBI-83353	1304-CP-7-25-L	2-5-2	4	2	Large	110
IBI-83354	1304-CP-7-25-XL	2-5-2	4	2	X-Large	110
IBI-83358	1304-CP-7-25-M-TH	2-5-2	4	2	Medium	110
IBI-83359	1304-CP-7-25-L-TH	2-5-2	4	2	Large	110
IBI-83360	1304-CP-7-25-XL-TH	2-5-2	4	2	X-Large	110
IBI-83362	1304-CP-7-25-L-TE2	2-5-2	2	2	Large	110
IBI-83376	1304-CP-7-25-FL	2-5-2	4	2	Flutter	110
IBI-83382	1304-CP-7-25-L1	2-5-2	4	2	Large1	110
IBI-83388	1304-CP-7-25-L1	2	4	2	Large1	110

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641



Therapy™ Cool Path™ Ablation Catheters

4 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- All braided shaft
- Double mesh stainless steel braiding to the distal tip
- Six equidistant ports on the distal tip
- Uni-directional steering
- Push/pull handle


Ordering Information

7 F irrigated ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
IBI-83502	1304-CP-7-25-L-AB	2-5-2	4	2	Large	110
IBI-83503	1304-CP-7-25-M-AB	2-5-2	4	2	Medium	110
IBI-83504	1304-CP-7-25-XL-AB	2-5-2	4	2	X-Large	110
IBI-83505	1304-CP-7-25-L1-AB	2-5-2	4	2	Large1	110
IBI-83506	1304-CP-7-25-FL-AB	2-5-2	4	2	Flutter	110

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641



Therapy™ Cool Path™ Duo

Ablation Catheters
4 mm Tip Thermocouple
Quadripolar
7 F

Product Highlights

- All braided shaft
- Double mesh stainless steel braiding to the distal tip
- Two rows of six equidistant ports offset by 30 degrees
- Proximal and distal tip electrode irrigation
- Uni-directional steering
- Push/pull handle


Ordering Information

7 F irrigated ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE
IBI-83562	1304-CP2-7-1.5(5)2-M-TE4BE1AB	1.5-5-2	4	1	Medium
IBI-83563	1304-CP2-7-1.5(5)2-L-TE4BE1AB	1.5-5-2	4	1	Large
IBI-83564	1304-CP2-7-1.5(5)2-XL-TE4BE1AB	1.5-5-2	4	1	X-Large
IBI-83565	1304-CP2-7-1.5(5)2-L1-TE4BE1AB	1.5-5-2	4	1	Large1

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641, IBI-85644



Therapy™ Cool Flex™

Ablation Catheter

Irrigated 4 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- Flexible tip electrode can conform to tissue
- Laser cut tip designed for optimal irrigation flow profile
- 0.5 mm spacing between tip and first band electrode to enhance signal quality
- Uni-directional steering
- Push/pull handle


Ordering Information

7 F irrigated ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
A088015	1304-CF-7-0.5(5)2-M-TE4BE1EB	0.5-5-2	4	1	Medium	110
A088016	1304-CF-7-0.5(5)2-L-TE4BE1EB	0.5-5-2	4	1	Large	110
A088017	1304-CF-7-0.5(5)2-XL-TE4BE1EB	0.5-5-2	4	1	X-Large	110
A088018	1304-CF-7-0.5(5)2-L1-TE4BE1EB	0.5-5-2	4	1	L1	110

Required Catheter Connecting Cable – Page A.36

SJM: IBI-85641



Contact™ Therapy™ Cool Path™ Duo

Ablation Catheter

Irrigated 4 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- Compatible with EnSite™ Contact™ Technology
- Uni-directional steering with push/pull handle
- Proximal and distal tip electrode irrigation
- Two rows of six equidistant ports offset by 30 degrees
- Double mesh stainless steel braiding to the distal tip


Ordering Information

7 F irrigated ablation catheter (1 unit per box)

REORDER NUMBER	DESCRIPTION	ELECTRODE SPACING (MM)	TIP ELECTRODE (MM)	BAND ELECTRODE (MM)	CURVE TYPE	USABLE LENGTH (CM)
A088019	1304-CP2-7-1.5(5)2-M-TE4BE1AB-CN	1.5-5-2	4	1	Medium	110
A088020	1304-CP2-7-1.5(5)2-L-TE4BE1AB-CN	1.5-5-2	4	1	Large	110
A088022	1304-CP2-7-1.5(5)2-L1-TE4BE1AB-CN	1.5-5-2	4	1	L1	110

Required Catheter Connecting Cables – Page A.36

A402892



Cool Point™

Irrigation Pump

Product Highlights

- Peristaltic pump
- Two-way communication with the IBI-1500T11 Cardiac Ablation Generator
- 2 microliter bubble detectors
- In-line pressure sensor to detect downstream occlusion
- Compact, lightweight design
- Ergonomic carrying handle
- Includes: Irrigation Pump, Power Cord, Pole Clamp, 1779 Connecting Cable, Operator's Manual

Ordering Information


REORDER NUMBER	DESCRIPTION
85784	Cool Point™ Irrigation Pump

Required Cool Point Tubing Set

REORDER NUMBER	DESCRIPTION	LENGTH (CM)
85785	Cool Point™ Tubing Set (sold individually)	329

Required Cool Point Accessories

REORDER NUMBER	DESCRIPTION	LENGTH (CM)
85786	Communications cable, 20-feet, between Cool Point™ Irrigation Pump and IBI-1500T11 Cardiac Ablation Generator	600



Ampere™ RF

Ablation Generator

Product Highlights

- Specifically designed to work alongside the EnSite™ Cardiac Mapping System, including an integrated GenConnect module
- Designed for compatibility with EP-WorkMate™ Recording System, WorkMate Claris™ Recording System and other major recording systems
- Choice of power or temperature control modes
- Large, colorful screen and buttons for easy reading and programming
- New programming options, including automatic irrigation flow adjustment and variable power ramp-up times
- Up to four presets for easy switching between settings
- Optional full-color screen remote control and footswitch available

Ordering Information

REORDER NUMBER	DESCRIPTION
H700495	Ampere™ Generator Kit (includes generator + footswitch, 2.5 m cable)
H700491	Ampere™ Remote Control (includes 15 m fiber cord)
H700498	Ampere™ Footswitch, 2.5 m
H700499	Ampere Footswitch, 5 m
H700500	Ampere Footswitch, 10 m
H700501	Ampere Footswitch, 15 m



Ampere™ RF Ablation Generator



Ampere™ Remote Control

TactiSys™ Quartz Equipment

Product Highlights

- Fully integrated display with the EnSite Precision™ Cardiac Mapping System and Contact Force Module
- Optical white light interferometry sensor based technology for enhanced resolution and measurement stability
- Fast sampling rate for accurate contact force information
- No lengthy startup calibration procedure required
- Proprietary indices of Force Time Integral (FTI™) and Lesion Index (LSI™)

Ordering Information

REORDER NUMBER	DESCRIPTION
PN-004 400	TactiSys™ Quartz Equipment
PN-004 515	TactiSys™ Quartz Equipment Ampere RF Cable
H702500	EnSite™ Contact Force Module v.2.0 for use with the EnSite Precision™ Cardiac Mapping System
TSM-ABV-GBL	TactiSys™ Quartz Equipment Mounting Bracket (for placement above EnSite Precision™ Link)
TSM-BSD-GBL	TactiSys™ Quartz Equipment Mounting Bracket (for placement beside EnSite Precision™ Link)
BED-CLP-WE	Bedrail Clamp (Western Europe)
BED-CLP-UK	Bedrail Clamp (UK)



Catheter Extension Cables

Livewire™ TC Ablation Catheters


Product Highlights

- Cable clearly identified for ablation
- Connection allows for temperature monitoring
- Guided channel for fast convenient connection


Ordering Information

Extension cable (1 unit per box)


REORDER NUMBER	LENGTH (CM)	RF GENERATOR TYPE
402505	300	Stockert Shuttle
402542	300	EPT-1000 XP (100 watt max)
402543	150	IBI-1500T Series (100 watt max)
402544	300	IBI-1500T Series (100 watt max)




IBI-1500T SERIES



EPT-1000 XP



ATAKR



STOCKERT SHUTTLE

Atakr is a trademark of Medtronic, Inc.

EPT-1000 XP is a trademark of Boston Scientific, Inc.

Catheter Connection Cables for

Therapy™ and Contact Enabled™
Ablation Catheters


Product Highlights

- Compatible cables for Therapy™ and FlexAbility™ catheters to compatible generators

Ordering Information

REORDER NUMBER	MODEL NUMBER	LENGTH (CM)	CATHETER	GENERATOR
IBI-85641	1641	250	Therapy™ Series Single Thermocouple/Therapy™ Dual-8™ / Therapy™ Cool Path™ / Therapy™ Cool Path™ Duo/ Therapy™ Cool Flex™ / FlexAbility™	Ampere™ RF Ablation Generator
IBI-85644	1641-A	250	Therapy™ Series Single Thermocouple/Therapy™ Dual-8™ / Therapy™ Cool Path™ / Therapy™ Cool Path™ Duo/ Therapy™ Cool Flex™	Ampere™ RF Ablation Generator
IBI-85708	1711-MA	250	Therapy™ Series Thermocouple	Medtronic Atakr
IBI-85709	1713-WA	250	Therapy™ Series Thermocouple	Stockert Shuttle
IBI-85711	1711-M	250	Therapy™ Series Thermocouple	Medtronic Atakr
IBI-85713	1713-W	250	Therapy™ Series, Single Thermocouple	Stockert Shuttle
IBI-85739	1739-W	250	Therapy™ Dual-8™	Stockert Shuttle
IBI-85761	1761-W	250	Therapy™ Series Thermistor; 4 mm tip	Stockert Shuffle
IBI-85763	1763-E	250	Therapy™ Series Thermistor; 4 mm tip	EPT-1000 XP
A402892	1611	250	Contact Enabled™	Ampere™ RF Ablation Generator


1711-M 1711-MA



1641, 1641-A, 1713-W, 1713-WA, 1739-W, 1761-W,
1763-E, 1763-EA



1611



Atakr is a trademark of Medtronic, Inc.
EPT-1000 TC is a trademark of Boston Scientific, Inc.

